

### **Fixation of Tissue Culture Cells for Fluorescence Microscopy:**

(unless otherwise specified, all solutions are prepared in PBS. **Preparation of solutions in water will result in osmotic lysis of cells).**

1. Seed cells at desired density on glass cover slides
2. Following day (or in 2 days), wash cells 1X with PBS
3. Fix for 15-20 min. in 4% paraformaldehyde'
  - wash 2X PBS
4. Quench with 50 mM NH<sub>4</sub>Cl (to get rid of background fluorescence from free aldehyde groups) 10 min.
  - Wash 2X PBS
5. Permeabilize 15' in 30 µg/mL digitonin
6. Wash 2X PBS
7. Block 30' RT (or O/N 4 °C) in 0.1% Teleost Fish gelatin/0.5% Triton X-100
  - Cells are ready to be treated with antibody

Rationale: digitonin is a gentle detergent, 0.5% triton X-100 will maintain permeabilization created by digitonin without generating new permeabilizations  
Gelatin vs. BSA – gelatin should bind more nonspecifically.