

Invasion/Replication Assay

(from Anthea Lee, Falkow Lab) modified by Igor Brodsky 11/12/98

Day 1:

1. Seed macrophages at 2.5×10^5 cells/well (24-well dishes)
Use 3 wells per data point (triplicate)
Hemocytometer: count cells in grid. #cells $\times 10^4$ = cells/ml of preparation used to count
2. Inoculate standing O/N (s. typhimurium) cultures into high salt LB

Day2:

1. Subculture O/N cultures (1:20-1:100) again standing – allow to grow up (optional – may give greater invasiveness)
Calculate cfu/ml from OD_{600} : OD_{600} of 1.0 = 1.2×10^9 cfu/ml
2. Aspirate medium from cells, wash 1X w/ PBS prewarmed to 37 °C. Aspirate PBS, Replace with 200 μ l of prewarmed fresh medium.
3. Infect cells with appropriate MOI (usually 10) in 10-50 μ l
Plate input (use as standard for % invasion)
4. Incubate abcteria with cells for 30-60' to allow invasion
5. Aspirate DMEM, wash 2x with DMEM, aspirate washes, add 1ml DMEM containing 100 μ g gentamycin /ml to each well.
6. Incubate with gentamycin for 90-120' to allow killing of extracellular bacteria
7. Wash 3x w/ DMEM

For Invasion Assay:

8. Aspirate medium, add 200 μ l 1% Triton X-100. Incubate RT 5' to allow host-cell lysis
9. Pipette vigorously up&down to ensure complete/uniform lysis
10. Add 800 ml DMEM
11. Make 10-fold serial dilutions, plate 10^{-3} and 10^{-4} (other dilutions if desired)

For Replication Assay:

8. Aspirate medium, add back DMEM containing 10 μ g/ml gentamycin to each well
9. At each desired time point, wash 1-2X with DMEM, follow steps 8-11 of *Invasion Assay*