

## Macrophage Invasion/phagocytosis Assay

Mφ: RAW264.7

Bacteria: Yersinia

Tissue culture media: DMEM + 10% DFCS + glutamine + 1mM Na Pyruvate

Assay:

1.  Mφ culture:  
RAW's seeded 12-24 h before experiment @  $2.5 \times 10^5$  mf/well OR  $1.5 \times 10^5$  for microscopy
2.  Bacterial cultures
  - A.  grown O/N @ 26°C shaking in 2xYT media
  - B.  O/N's subcultured 1:50 in 2.5 ml of 2xYT + 20mM NaOxalate + 20mM MgCl<sub>2</sub>  
  
per 10ml = 9 ml 2xYT  
0.8 ml NaOxalate (200mM)  
0.2 ml MgCl<sub>2</sub> (1M)
  - C.  shake @ 26°C for 2 h followed by a shift to 37°C shaking for 2 h

Invasion

- A.  add approx.  $2.5 \times 10^6$  CFU/well of mφ (read OD<sub>600</sub> 1.0 =  $2 \times 10^9$  CFU/ml)
- B.  centrifuge bacteria onto monolayer @ 2000 rpm for 5 min @ RT
- C.  place infected monolayers @ 37°C + 5% CO<sub>2</sub> for 90 min OR do time course and fix accordingly.
- D.  Aspirate media and gently add 1 ml of pre-warmed DMEM containing 100 µg/ml gentamicin to each sample.
- E.  Return to incubator for 90 min to kill extracellular bacteria
- F.  Aspirate gentamicin containing media and wash wells 1X with 0.5 ml pre-warmed DMEM
- G.  Aspirate and add 200 µl of 1% triton X-100 to each well and allow mφ to lyse for 5 min @ RT
- H.  Add 800 µl of LB media to each well and titer sample wells appropriately to determine output CFU